

For Research Use Only

# TFE3 Monoclonal antibody, PBS Only

Catalog Number: 67319-1-PBS



## Basic Information

Catalog Number:

67319-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6320

GenBank Accession Number:

BC026027

GeneID (NCBI):

7030

UNIPROT ID:

P19532

Full Name:

transcription factor binding to IGHM  
enhancer 3

Calculated MW:

62 kDa

Observed MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

1E3D7

## Applications

Tested Applications:

WB, Indirect ELISA, IHC

Species Specificity:

Human, mouse, rat, pig

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67319-1-Ig (TFE3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67319-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human malignant melanoma tissue slide using 67319-1-Ig (TFE3 antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 67319-1-PBS in a different storage buffer formulation.