

For Research Use Only

HSPB6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67327-2-PBS



Basic Information

Catalog Number:

67327-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10484

GenBank Accession Number:

BC068046

GeneID (NCBI):

126393

UNIPROT ID:

O14558

Full Name:

heat shock protein, alpha-crystallin-related, B6

Calculated MW:

160 aa, 17 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A6G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50552-1, HSPB6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67327-2-PBS. Detection antibody: 67327-3-PBS. Standard: Ag10484. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50552-4, HSPB6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67327-2-PBS. Detection antibody: 67327-6-PBS. Standard: Ag10484. Range: 6.25-100 ng/mL.