

For Research Use Only

# NICN1 Monoclonal antibody, PBS Only



Catalog Number: 67328-1-PBS

## Basic Information

Catalog Number: 67328-1-PBS	GenBank Accession Number: BC050653	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84276	CloneNo.: 5B7G6
Source: Mouse	UNIPROT ID: Q9B5H3	
Isotype: IgG2b	Full Name: nicolin 1	
Immunogen Catalog Number: AG6581	Calculated MW: 24 kDa	
	Observed MW: 27 kDa	

## Applications

Tested Applications:  
WB, Indirect ELISA

Species Specificity:  
Human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



U-251 cells were subjected to SDS PAGE followed by western blot with 67328-1-Ig (NICN1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67328-1-PBS in a different storage buffer formulation.