

For Research Use Only

ACSM5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67334-2-PBS



Basic Information

Catalog Number:

67334-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10056

GenBank Accession Number:

BC016703

GeneID (NCBI):

54988

UNIPROT ID:

Q6NUN0

Full Name:

acyl-CoA synthetase medium-chain
family member 5

Calculated MW:

208 aa, 23 kDa

Purification Method:

Protein G purification

CloneNo.:

2A10B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50791-1, ACSM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67334-2-PBS. Detection antibody: 67334-3-PBS. Standard: Ag9373. Range: 1.563-100 ng/mL.