

For Research Use Only

FUT3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67344-2-PBS



Basic Information

Catalog Number: 67344-2-PBS	GenBank Accession Number: BC108675	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2525	CloneNo.: 2E2B4
Source: Mouse	UNIPROT ID: P21217	
Isotype: IgG1	Full Name: fucosyltransferase 3 (galactoside 3(4)-L-fucosyltransferase, Lewis blood group)	
Immunogen Catalog Number: AG12581	Calculated MW: 361 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50464-1, FUT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67344-2-PBS. Detection antibody: 67344-3-PBS. Standard: Ag12581. Range: 0.098-12.5 ng/mL