

For Research Use Only

PTGS1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67346-5-PBS



Basic Information

Catalog Number:

67346-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28280

GenBank Accession Number:

BC029840

GeneID (NCBI):

5742

UNIPROT ID:

P23219

Full Name:

prostaglandin-endoperoxide synthase

1 (prostaglandin G/H synthase and

cyclooxygenase)

Calculated MW:

599 aa, 69 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E8E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50848-2, PTGS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67346-4-PBS. Detection antibody: 67346-5-PBS. Standard:Ag28280. Range: 0.195-50 ng/mL



Cytometric bead array standard curve of MP50848-3, PTGS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67346-5-PBS. Detection antibody: 67346-3-PBS. Standard:Ag28280. Range: 0.391-100 ng/mL