

For Research Use Only

# SEC63 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67352-2-PBS



## Basic Information

Catalog Number:

67352-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5783

GenBank Accession Number:

BC047221

GeneID (NCBI):

11231

UNIPROT ID:

Q9UGP8

Full Name:

SEC63 homolog (S. cerevisiae)

Calculated MW:

88 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A3B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50799-1, SEC63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67352-2-PBS. Detection antibody: 67352-3-PBS. Standard: Ag5783. Range: 0.195-12.5 ng/mL