

For Research Use Only

RALBP1 Monoclonal antibody, PBS Only

Catalog Number: 67359-1-PBS



Basic Information

Catalog Number: 67359-1-PBS	GenBank Accession Number: BC013126	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 10928	CloneNo.: 1B2F5
Source: Mouse	UNIPROT ID: Q15311	
Isotype: IgG1	Full Name: raIA binding protein 1	
Immunogen Catalog Number: AG28618	Calculated MW: 655 aa, 76 kDa	
	Observed MW: 95 kDa	

Applications

Tested Applications:
WB, FC (Intra), Indirect ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 67359-1-Ig (RALBP1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67359-1-PBS in a different storage buffer formulation.

1x10⁶ HEK-293 cells were intracellularly stained with 0.25 ug RALBP1 Monoclonal antibody (67359-1-Ig, Clone:1B2F5) and CoraLite®488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.25 ug Mouse IgG1 isotype control Mouse McAb (66360-1-Ig, Clone: 1F8D3) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 67359-1-PBS in a different storage buffer