

For Research Use Only

FBXO11 Monoclonal antibody, PBS Only



Catalog Number: 67365-1-PBS

Basic Information

Catalog Number: 67365-1-PBS	GenBank Accession Number: BC012728	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 80204	CloneNo.: 3C5E1
Source: Mouse	UNIPROT ID: Q86XK2	
Isotype: IgG1	Full Name: F-box protein 11	
Immunogen Catalog Number: AG24401	Calculated MW: 927 aa, 104 kDa	
	Observed MW: 94-104 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67365-1-Ig (FBXO11 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67365-1-PBS in a different storage buffer formulation.