

For Research Use Only

CAPN3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67366-3-PBS



Basic Information

Catalog Number:

67366-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13179

GenBank Accession Number:

BC004883

GeneID (NCBI):

825

UNIPROT ID:

P20807

Full Name:

calpain 3, (p94)

Calculated MW:

18 kDa

Purification Method:

Protein G purification

CloneNo.:

2F3B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51120-1, CAPN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67366-2-PBS. Detection antibody: 67366-3-PBS. Standard:Ag13179. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51120-2, CAPN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67366-4-PBS. Detection antibody: 67366-3-PBS. Standard:Ag13179. Range: 0.781-100 ng/mL