

For Research Use Only

CEP164 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67370-2-PBS



Basic Information

Catalog Number: 67370-2-PBS	GenBank Accession Number: BC000602	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 22897	CloneNo.: 3B7F2
Source: Mouse	UNIPROT ID: Q9UPV0	
Isotype: IgG1	Full Name: centrosomal protein 164kDa	
Immunogen Catalog Number: AG17699	Calculated MW: 1460 aa, 164 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50474-1, CEP164 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67370-2-PBS. Detection antibody: 67370-3-PBS. Standard: Ag17699. Range: 0.195-50 ng/mL



Cytometric bead array standard curve of MP50474-2, CEP164 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67370-2-PBS. Detection antibody: 67370-4-PBS. Standard: Ag17699. Range: 0.098-50 ng/mL