

For Research Use Only

RBMS2 Monoclonal antibody

Catalog Number: 67395-1-Ig



Basic Information

Catalog Number:

67395-1-Ig

Size:

1800 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27721

GenBank Accession Number:

BC027863

GeneID (NCBI):

5939

UNIPROT ID:

Q15434

Full Name:

RNA binding motif, single stranded interacting protein 2

Calculated MW:

407 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

1E5C5

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB: NCI-H1299 cells, U2OS cells, LNCaP cells, A549 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67395-1-Ig (RBMS2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.