

For Research Use Only

# GATAD2B Monoclonal antibody

Catalog Number: 67396-1-Ig



## Basic Information

Catalog Number: 67396-1-Ig	GenBank Accession Number: BC069419	Purification Method: Protein G purification
Size: Source: Mouse	GeneID (NCBI): 57459	CloneNo.: 2C3F5
Isotype: IgG1	UNIPROT ID: Q8WXI9	
Immunogen Catalog Number: AG22615	Full Name: GATA zinc finger domain containing 2B	

## Applications

Tested Applications:  
ELISA

Species Specificity:  
Human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C.

Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data