

For Research Use Only

# SEMA7A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67397-3-PBS

## Basic Information

Catalog Number:

67397-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28804

GenBank Accession Number:

BC101643

GeneID (NCBI):

8482

UNIPROT ID:

O75326

Full Name:

semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Calculated MW:

666 aa, 75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F5B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50321-1, SEMA7A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67397-2-PBS. Detection antibody: 67397-3-PBS. Standard: Ag28804. Range: 0.098-25 ng/mL