

For Research Use Only

SEMA7A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67397-5-PBS

Basic Information

Catalog Number:

67397-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28804

GenBank Accession Number:

BC101643

GeneID (NCBI):

8482

UNIPROT ID:

O75326

Full Name:

semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Calculated MW:

666 aa, 75 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A9D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50321-2, SEMA7A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67397-4-PBS. Detection antibody: 67397-5-PBS. Standard:SEMA7A. Range: 0.391-25 ng/mL