

For Research Use Only

TYK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67411-2-PBS



Basic Information

| | | |
|--|--|--|
| Catalog Number: 67411-2-PBS | GenBank Accession Number: BC014243 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 7297 | CloneNo.: 1H5B4 |
| Source: Mouse | UNIPROT ID: P29597 | |
| Isotype: IgG1 | Full Name: tyrosine kinase 2 | |
| Immunogen Catalog Number: AG9442 | Calculated MW: 1187 aa, 134 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51012-1, TYK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67411-2-PBS. Detection antibody: 67411-3-PBS. Standard: Ag9442. Range: 0.098-100 ng/mL.