

For Research Use Only

TYK2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67411-3-PBS



Basic Information

Catalog Number: 67411-3-PBS	GenBank Accession Number: BC014243	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 7297	CloneNo.: 1C10F4
Source: Mouse	UNIPROT ID: P29597	
Isotype: IgG2a	Full Name: tyrosine kinase 2	
Immunogen Catalog Number: AG9442	Calculated MW: 1187 aa, 134 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51012-1, TYK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67411-2-PBS. Detection antibody: 67411-3-PBS. Standard: Ag9442. Range: 0.098-100 ng/mL.