

For Research Use Only

TYK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67411-4-PBS



Basic Information

Catalog Number:

67411-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9442

GenBank Accession Number:

BC014243

GeneID (NCBI):

7297

UNIPROT ID:

P29597

Full Name:

tyrosine kinase 2

Calculated MW:

1187 aa, 134 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H1G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51012-2, TYK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67411-4-PBS. Detection antibody: 67411-5-PBS. Standard: Ag9442. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51012-3, TYK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67411-4-PBS. Detection antibody: 67411-6-PBS. Standard: Ag9442. Range: 0.391-100 ng/mL.