

For Research Use Only

ERO1L Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67416-4-PBS



Basic Information

Catalog Number: 67416-4-PBS	GenBank Accession Number: BC008674	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 30001	CloneNo.: 2B5E5
Source: Mouse	UNIPROT ID: Q96HE7	
Isotype: IgG1	Full Name: ERO1-like (S. cerevisiae)	
Immunogen Catalog Number: AG29910	Calculated MW: 468 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51017-2, ERO1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67416-2-PBS. Detection antibody: 67416-4-PBS. Standard: Ag29910. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP51017-3, ERO1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67416-4-PBS. Detection antibody: 67416-5-PBS. Standard: Ag29910. Range: 0.781-50 ng/mL.