

For Research Use Only

CYP2W1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67419-4-PBS



Basic Information

Catalog Number:

67419-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29982

GenBank Accession Number:

BC025761

GeneID (NCBI):

54905

UNIPROT ID:

Q8TAV3

Full Name:

cytochrome P450, family 2, subfamily W, polypeptide 1

Calculated MW:

54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C1G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50810-2, CYP2W1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67419-2-PBS. Detection antibody: 67419-4-PBS. Standard: Ag29982. Range: 0.781-200 ng/mL.