

For Research Use Only

SAMM50 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67425-2-PBS



Basic Information

Catalog Number:

67425-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15155

GenBank Accession Number:

BC007830

GeneID (NCBI):

25813

UNIPROT ID:

Q9Y512

Full Name:

sorting and assembly machinery component 50 homolog (S. cerevisiae)

Calculated MW:

469 aa, 52 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B1G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50503-1, SAMM50 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67425-2-PBS. Detection antibody: 67425-3-PBS. Standard: Ag15155. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50503-3, SAMM50 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67425-4-PBS. Detection antibody: 67425-2-PBS. Standard: Ag15155. Range: 0.098-100 ng/mL.