

For Research Use Only

# SAMM50 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67425-4-PBS



## Basic Information

Catalog Number:

67425-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG15155

GenBank Accession Number:

BC007830

GeneID (NCBI):

25813

UNIPROT ID:

Q9Y512

Full Name:

sorting and assembly machinery component 50 homolog (S. cerevisiae)

Calculated MW:

469 aa, 52 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A8E2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50503-2, SAMM50 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67425-4-PBS. Detection antibody: 67425-3-PBS. Standard:Ag15155. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50503-3, SAMM50 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67425-4-PBS. Detection antibody: 67425-2-PBS. Standard:Ag15155. Range: 0.098-100 ng/mL