

For Research Use Only

PDPN Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67432-3-PBS

Basic Information

Catalog Number: 67432-3-PBS	GenBank Accession Number: BC022812	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10630	CloneNo.: 1B12E4
Source: Mouse	UNIPROT ID: Q86YL7	
Isotype: IgG2a	Full Name: podoplanin	
Immunogen Catalog Number: AG17691	Calculated MW: 238 aa, 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50340-2, PDPN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67432-3-PBS. Detection antibody: 67432-4-PBS. Standard: Ag17691. Range: 3.125-100 ng/mL.