

For Research Use Only

LIPG Monoclonal antibody, PBS Only

Catalog Number: 67434-1-PBS



Basic Information

Catalog Number: 67434-1-PBS	GenBank Accession Number: BC060825	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9388	CloneNo.: 1B10B10
Source: Mouse	UNIPROT ID: Q9Y5X9	
Isotype: IgG1	Full Name: lipase, endothelial	
Immunogen Catalog Number: AG29482	Calculated MW: 57 kDa	
	Observed MW: 57 and 42 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, IHC

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67434-1-ig (LIPG antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67434-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 67434-1-ig (LIPG antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67434-1-PBS in a different storage buffer formulation.