

For Research Use Only

GPR177 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67435-5-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 67435-5-PBS | GenBank Accession Number: BC110826 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 79971 | CloneNo.: 1D1H8 |
| Source: Mouse | UNIPROT ID: Q5T9L3 | |
| Isotype: IgG1 | Full Name: G protein-coupled receptor 177 | |
| Immunogen Catalog Number: AG12446 | Calculated MW: 541 aa, 62 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50915-2, GPR177 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67435-4-PBS. Detection antibody: 67435-5-PBS. Standard: Ag12446. Range: 6.25-200 ng/mL