

For Research Use Only

RAGE Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67437-1-PBS



Basic Information

Catalog Number: 67437-1-PBS	GenBank Accession Number: BC053536	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 5891	CloneNo.: 2C8D6
Source: Mouse	UNIPROT ID: Q9UQ07	
Isotype: IgG1	Full Name: renal tumor antigen	
Immunogen Catalog Number: AG29406	Calculated MW: 231 aa, 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50901-1, RAGE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67437-1-PBS. Detection antibody: 67437-2-PBS. Standard:Ag29406. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50901-2, RAGE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67437-1-PBS. Detection antibody: 67437-3-PBS. Standard:Ag29406. Range: 0.781-100 ng/mL