

For Research Use Only

SOX9 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67439-2-PBS



Basic Information

Catalog Number:

67439-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28620

GenBank Accession Number:

BC007951

GeneID (NCBI):

6662

UNIPROT ID:

P48436

Full Name:

SRV (sex determining region Y)-box 9

Calculated MW:

509 aa, 56 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H3C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51117-1, SOX9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67439-2-PBS. Detection antibody: 67439-3-PBS. Standard:Ag28620. Range: 0.781-50 ng/mL