

For Research Use Only

SOX9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67439-3-PBS



Basic Information

Catalog Number: 67439-3-PBS	GenBank Accession Number: BC007951	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6662	CloneNo.: 3F7A3
Source: Mouse	UNIPROT ID: P48436	
Isotype: IgG1	Full Name: SRV (sex determining region Y)-box 9	
Immunogen Catalog Number: AG28620	Calculated MW: 509 aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51117-1, SOX9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67439-2-PBS. Detection antibody: 67439-3-PBS. Standard:Ag28620. Range: 0.781-50 ng/mL