

For Research Use Only

# NUDT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67443-4-PBS



## Basic Information

Catalog Number:

67443-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9802

GenBank Accession Number:

BC014618

GeneID (NCBI):

4521

UNIPROT ID:

P36639

Full Name:

nudix (nucleoside diphosphate linked moiety X)-type motif 1

Calculated MW:

179aa, 20 kDa; 197aa, 23 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C9F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50783-2, NUDT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67443-1-PBS. Detection antibody: 67443-4-PBS. Standard: Ag9802. Range: 0.195-100 ng/mL.