

For Research Use Only

CD62L Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67444-2-PBS



Basic Information

Catalog Number: 67444-2-PBS	GenBank Accession Number: BC020758	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6402	CloneNo.: 3C3B8
Source: Mouse	ENSEMBL Gene ID: ENSG00000188404	
Isotype: IgG1	UNIPROT ID: P14151	
Immunogen Catalog Number: AG29297	Full Name: selectin L	
	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50444-1, CD62L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67444-2-PBS. Detection antibody: 67444-3-PBS. Standard: Ag29297. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50444-2, CD62L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67444-2-PBS. Detection antibody: 67444-4-PBS. Standard: Ag29297. Range: 0.098-6.25 ng/mL