

For Research Use Only

CD62L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67444-4-PBS



Basic Information

Catalog Number:

67444-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29297

GenBank Accession Number:

BC020758

GeneID (NCBI):

6402

ENSEMBL Gene ID:

ENSG00000188404

UNIPROT ID:

P14151

Full Name:

selectin L

Calculated MW:

42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A4B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50444-2, CD62L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67444-2-PBS. Detection antibody: 67444-4-PBS. Standard:Ag29297. Range: 0.098-6.25 ng/mL