

For Research Use Only

MYC Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67447-3-PBS

Basic Information

Catalog Number: 67447-3-PBS	GenBank Accession Number: BC000141	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4609	CloneNo.: 3D2A7
Source: Mouse	UNIPROT ID: P01106	
Isotype: IgG1	Full Name: v-myc myelocytomatosis viral oncogene homolog (avian)	
Immunogen Catalog Number: AG17519	Calculated MW: 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50347-1, MYC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67447-2-PBS. Detection antibody: 67447-3-PBS. Standard:Ag17519. Range: 0.781-50 ng/mL