

For Research Use Only

OMA1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67449-2-PBS



Basic Information

Catalog Number: 67449-2-PBS	GenBank Accession Number: BC012915	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 115209	CloneNo.: 2E1D12
Source: Mouse	UNIPROT ID: Q96E52	
Isotype: IgG1	Full Name: OMA1 homolog, zinc metallopeptidase (S. cerevisiae)	
Immunogen Catalog Number: AG11201	Calculated MW: 524 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50571-2, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-3-PBS. Detection antibody: 67449-2-PBS. Standard: Ag11201. Range: 0.781-100 ng/mL.