

For Research Use Only

ME2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67457-4-PBS



Basic Information

Catalog Number:

67457-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG20741

GenBank Accession Number:

BC000147

GeneID (NCBI):

4200

UNIPROT ID:

P23368

Full Name:

malic enzyme 2, NAD(+)-dependent,
mitochondrial

Calculated MW:

65 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

4F9C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50528-1, ME2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67457-2-PBS. Detection antibody: 67457-4-PBS. Standard:Ag20741. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50528-2, ME2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67457-3-PBS. Detection antibody: 67457-4-PBS. Standard:Ag20741. Range: 0.781-25 ng/mL.