

For Research Use Only

ZNF746 Monoclonal antibody, PBS Only



Catalog Number: 67463-1-PBS

Basic Information

Catalog Number: 67463-1-PBS	GenBank Accession Number: BC068505	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 155061	CloneNo.: 3A12C5
Source: Mouse	UNIPROT ID: Q6NUN9	
Isotype: IgG1	Full Name: zinc finger protein 746	
Immunogen Catalog Number: AG21582	Calculated MW: 644 aa, 69 kDa	
	Observed MW: 85-90 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

ZNF746, also named as PARIS, which belongs to the krueppel C2H2-type zinc-finger protein family. ZNF746 mainly localizes to the cytoplasm and probably translocates to the nucleus to repress selected genes. ZNF746 as a transcription repressor that specifically binds to the 5'-TATTTT[T/G]-3' consensus sequence on promoters and repress transcription of PGC-1-alpha (PPARGC1A), thereby playing a role in regulation of neuron death. ZNF746 may act as a downstream effector of PRKN and contribute to neurodegeneration in Parkinson disease cases caused by defects in PRKN: its accumulation due to the absence of PRKN, followed by up-regulation of PPARGC1A, could lead to the selective loss of dopamine neurons in the substantia nigra (PMID: 28122242). The calculated molecular weight of ZNF746 is 69 kDa, but Phosphorylation of ZNF746 is about 80-85 kDa (PMID: 28122242).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 67463-1-Ig (ZNF746 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67463-1-PBS in a different storage buffer formulation.