

For Research Use Only

MGME1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67468-2-PBS



Basic Information

Catalog Number:

67468-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19644

GenBank Accession Number:

BC016869

GeneID (NCBI):

92667

UNIPROT ID:

Q9BQP7

Full Name:

chromosome 20 open reading frame

72

Calculated MW:

344 aa, 39 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F7A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50593-1, MGME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67468-2-PBS. Detection antibody: 67468-3-PBS. Standard: Ag19644. Range: 12.5-200 ng/mL



Cytometric bead array standard curve of MP50593-2, MGME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67468-2-PBS. Detection antibody: 67468-4-PBS. Standard: Ag19644. Range: 6.25-200 ng/mL