

For Research Use Only

MGME1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67468-3-PBS



Basic Information

Catalog Number: 67468-3-PBS	GenBank Accession Number: BC016869	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 92667	CloneNo.: 1H3D1
Source: Mouse	UNIPROT ID: Q9BQP7	
Isotype: IgG1	Full Name: chromosome 20 open reading frame 72	
Immunogen Catalog Number: AG19644	Calculated MW: 344 aa, 39 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50593-1, MGME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67468-2-PBS. Detection antibody: 67468-3-PBS. Standard: Ag19644. Range: 12.5-200 ng/mL