

For Research Use Only

# PIP5K1A Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67483-2-PBS



## Basic Information

Catalog Number:

67483-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8355

GenBank Accession Number:

BC007833

GeneID (NCBI):

8394

UNIPROT ID:

Q99755

Full Name:

phosphatidylinositol-4-phosphate 5-kinase, type I, alpha

Calculated MW:

56 kDa, 63 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E9C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50903-1, PIP5K1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67483-2-PBS. Detection antibody: 67483-3-PBS. Standard: Ag8355. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50903-3, PIP5K1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67483-6-PBS. Detection antibody: 67483-2-PBS. Standard: Ag8355. Range: 0.195-100 ng/mL.