

For Research Use Only

# AHCYL2 Monoclonal antibody, PBS Only



Catalog Number: 67484-1-PBS

## Basic Information

<b>Catalog Number:</b> 67484-1-PBS	<b>GenBank Accession Number:</b> BC024325	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 23382	<b>CloneNo.:</b> 1D1F6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q96HN2	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> S-adenosylhomocysteine hydrolase-like 2	
<b>Immunogen Catalog Number:</b> AG19614	<b>Calculated MW:</b> 611 aa, 67 kDa	
	<b>Observed MW:</b> 67 and 57 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA, IHC

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 67484-1-Ig (AHCYL2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67484-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 67484-1-Ig (AHCYL2 antibody) at dilution of 1:600 (under 10x lens).. This data was developed using the same antibody clone with 67484-1-PBS in a different storage buffer formulation.