

For Research Use Only

S100P Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67485-4-PBS



Basic Information

Catalog Number: 67485-4-PBS	GenBank Accession Number: BC006819	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6286	CloneNo.: 1F7G1
Source: Mouse	UNIPROT ID: P25815	
Isotype: IgG1	Full Name: S100 calcium binding protein P	
Immunogen Catalog Number: AG19115	Calculated MW: 95 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51168-2, S100P Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67485-1-PBS. Detection antibody: 67485-4-PBS. Standard:Ag19115. Range: 6.25-100 ng/mL