

For Research Use Only

SPATA19 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67486-1-PBS



Basic Information

Catalog Number: 67486-1-PBS	GenBank Accession Number: BC058039	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 219938	CloneNo.: 1G12A11
Source: Mouse	UNIPROT ID: Q7Z5L4	
Isotype: IgG1	Full Name: spermatogenesis associated 19	
Immunogen Catalog Number: AG9808	Calculated MW: 167 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50970-1, SPATA19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67486-1-PBS. Detection antibody: 67486-2-PBS. Standard:Ag9808. Range: 6.25-200 ng/mL