

For Research Use Only

MFN2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67487-4-PBS



Basic Information

Catalog Number: 67487-4-PBS	GenBank Accession Number: BC017061	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9927	CloneNo.: 1D2B1
Source: Mouse	ENSEMBL Gene ID: ENSG00000116688	
Isotype: IgG1	UNIPROT ID: O95140	
Immunogen Catalog Number: AG29873	Full Name: mitofusin 2	
	Calculated MW: 757 aa, 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50442-2, MFN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67487-4-PBS. Detection antibody: 67487-5-PBS. Standard:Ag29873. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50442-3, MFN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67487-4-PBS. Detection antibody: 67487-6-PBS. Standard:Ag29873. Range: 1.563-100 ng/mL