

For Research Use Only

HRD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67488-2-PBS



Basic Information

Catalog Number: 67488-2-PBS	GenBank Accession Number: BC030530	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 84447	CloneNo.: 1D8C10
Source: Mouse	UNIPROT ID: Q86TM6	
Isotype: IgG1	Full Name: synovial apoptosis inhibitor 1, synoviolin	
Immunogen Catalog Number: AG4456	Calculated MW: 617 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50924-1, HRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67488-2-PBS. Detection antibody: 67488-3-PBS. Standard:Ag4456. Range: 0.391-100 ng/mL.