

For Research Use Only

# HRD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67488-3-PBS



## Basic Information

Catalog Number:

67488-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG4456

GenBank Accession Number:

BC030530

GeneID (NCBI):

84447

UNIPROT ID:

Q86TM6

Full Name:

synovial apoptosis inhibitor 1,  
synoviolin

Calculated MW:

617 aa, 68 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2G7B11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50924-1, HRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67488-2-PBS. Detection antibody: 67488-3-PBS. Standard: Ag4456. Range: 0.391-100 ng/mL.