

For Research Use Only

# TRIM5 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67492-2-PBS

## Basic Information

Catalog Number: 67492-2-PBS	GenBank Accession Number: BC021258	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 85363	CloneNo.: 2D12A10
Source: Mouse	UNIPROT ID: Q9C035	
Isotype: IgG1	Full Name: tripartite motif-containing 5	
Immunogen Catalog Number: AG30059	Calculated MW: 493 aa, 56 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50364-1, TRIM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67492-2-PBS. Detection antibody: 67492-3-PBS. Standard: Ag30059. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50364-2, TRIM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67492-2-PBS. Detection antibody: 67492-4-PBS. Standard: Ag30059. Range: 0.098-50 ng/mL.