

For Research Use Only

TRIM5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67492-4-PBS

Basic Information

Catalog Number: 67492-4-PBS	GenBank Accession Number: BC021258	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 85363	CloneNo.: 2E12E8
Source: Mouse	UNIPROT ID: Q9C035	
Isotype: IgG1	Full Name: tripartite motif-containing 5	
Immunogen Catalog Number: AG30059	Calculated MW: 493 aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50364-2, TRIM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67492-2-PBS. Detection antibody: 67492-4-PBS. Standard: Ag30059. Range: 0.098-50 ng/mL