

For Research Use Only

RAB2A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67501-3-PBS



Basic Information

Catalog Number:

67501-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7939

GenBank Accession Number:

BC008929

GeneID (NCBI):

5862

UNIPROT ID:

P61019

Full Name:

RAB2A, member RAS oncogene family

Calculated MW:

24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E3E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP51178-1, RAB2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67501-2-PBS. Detection antibody: 67501-3-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP51178-1, RAB2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67501-2-PBS. Detection antibody: 67501-3-PBS. Standard:Ag7939. Range: 3.125-100 ng/mL.