

For Research Use Only

RHOBTB3 Monoclonal antibody, PBS Only

Catalog Number: 67502-1-PBS



Basic Information

Catalog Number: 67502-1-PBS	GenBank Accession Number: BC041337	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 22836	CloneNo.: 1H2E3
Source: Mouse	UNIPROT ID: O94955	
Isotype: IgG1	Full Name: Rho-related BTB domain containing 3	
Immunogen Catalog Number: AG5782	Calculated MW: 611 aa, 69 kDa	
	Observed MW: 69 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 67502-1-Ig (RHOBTB3 antibody) at dilution of 1:300 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67502-1-PBS in a different storage buffer formulation.



K-562 cells were subjected to SDS PAGE followed by western blot with 67502-1-Ig (RHOBTB3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67502-1-PBS in a different storage buffer formulation.