

For Research Use Only

FES Monoclonal antibody

Catalog Number: 67524-1-Ig



Basic Information

Catalog Number: 67524-1-Ig	GenBank Accession Number: BC035357	Purification Method: Protein A purification
Size: 1600 µg/ml	GeneID (NCBI): 2242	CloneNo.: 1E1B11
Source: Mouse	UNIPROT ID: P07332	Recommended Dilutions: WB 1:2000-1:10000 IF 1:150-1:600
Isotype: IgG2b	Full Name: feline sarcoma oncogene	
Immunogen Catalog Number: AG28108	Calculated MW: 822 aa, 94 kDa	
	Observed MW: 93 kDa, 79 kDa, 87 kDa, 85 kDa	

Applications

Tested Applications: IF/ICC, WB, ELISA	Positive Controls:
Species Specificity: Human	WB : THP-1 cells, TF-1 cells, HL-60 cells, U-937 cells IF : A549 cells,

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67524-1-Ig (FES antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using 67524-1-Ig (FES antibody), at dilution of 1:300 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).