

For Research Use Only

# FES Monoclonal antibody, PBS Only

Catalog Number: 67524-1-PBS



## Basic Information

Catalog Number:

67524-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28108

GenBank Accession Number:

BC035357

GeneID (NCBI):

2242

UNIPROT ID:

P07332

Full Name:

feline sarcoma oncogene

Calculated MW:

822 aa, 94 kDa

Observed MW:

93 kDa, 79 kDa, 87 kDa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

1E1B11

## Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67524-1-Ig (FES antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67524-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using 67524-1-Ig (FES antibody), at dilution of 1:300 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67524-1-PBS in a different storage buffer formulation.